
# Test Cell ElyFlow

Elektrochemische Messzelle PTFE  
Electrochemical Test Cell PTFE

Artikelnummer / Item number 83400

User Manual

Version 2.0

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Lilienthalstrasse 146

D-34123 Kassel

+49 (0) 561 59190

info@gaskatel.de

[www.gaskatel.com](http://www.gaskatel.com)

## Languages

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# Operating instructions

## Table of contents

Overview .....	2
Scope of delivery .....	2
Declaration of conformity.....	3
General information.....	3
About these instructions.....	3
Intended use .....	3
Safety instructions.....	4
Use .....	5
Device overview .....	5
Operation.....	15
Operating conditions .....	23
Technical data.....	23
Disposal .....	24
Exclusion of liability .....	24

## Overview

### Scope of delivery

- 1 ElyFlow PTFE with
  - 1 x Counter electrode holder PTFE heatable (83410) with
    - 1 x Heating element (80258) with
      - 2 x Contact sockets 4 mm (80183)
    - 1 x Counter electrode Nickel for ElyFlow (80178)
    - 2 x Multilam plug M4 4 mm (89029)
    - 3 x Pin (80029)
  - 1 Working electrode holder suitable for gas supply (83420)
    - 2 x Straight connection ECTFE G1/8" (80013)
    - 2 x Multilam plug (89029)
    - 1 O-Ring EPDM 39.69 x 3.53 mm (80560)
    - 1 x Screw plug PTFE G1/8" (80031)
  - 1 x Analyte compartment PTFE with
    - 2 x Screw plug PTFE G1/8" (80031)
    - 4 x Screw plug PP G1/8" black (80133)
    - 1 x O-Ring EPDM 53,57x3,53 mm (80066)
    - 1 x O-Ring EPDM 39,69x3,53 mm (80560)
    - 6 x Nafion tube (83010)
    - 3 x Pin (80029)
    - 4 x Thumb screw M8x80 PA 6.6 (83005)
  - 1 x Counter electrode MMO (80182)
  - 4 x Thumb screw M8x50 PA6.6 (83007)
  - 2 x Test lead red 4 mm 1 mm<sup>2</sup> (89017)
  - 2 x Test lead black 4 mm 1 mm<sup>2</sup> (89019)
  - 4 x Connection coupling red 4 mm (89037)
  - 1 x Connection coupling black 4 mm (89038)
  - 1 x Hexagon socket spanner 6 x125mm 6kt. (83011)
  - 1 x Manual

## Optionally available accessories

Hydrogen reference electrode Mini-HydroFlex 81020

Analyte compartment PTFE 83430

Working electrode holder PMMA 83420-PMMA

Control box TP 83700

Temperature control box 83900

Electrolyte pump with power supply 83721

## Declaration of conformity

 This product complies with all EU directives applicable to this product.

## General information

### About these instructions

Read these instructions carefully. To ensure a long service life and reliable use, all instructions mentioned in this manual must be observed and followed.

Please keep these instructions at hand.

You can access the instructions and further information at [www.gaskatel.com](http://www.gaskatel.com).

### Intended use

Please consider the Good Laboratory Practice (GLP).

The test cell ElyFlow is suitable for electrochemical measurements on gas diffusion electrodes and other flat samples such as metals.

The test cells are designed for the typical three-electrode set-up. A counter electrode is already integrated.

The test cell ElyFlow can be heated via an integrated heating element. They may only be operated with the **Gaskatel temperature control box 83900** or **Gaskatel Control Box TP 83700**.

The test cell ElyFlow may be operated up to **100°C**. In case of using the working electrode holder made of PMMA, the cell may be operated up to 80°C.

The **filling level of the electrolyte must be checked regularly**, as the level can change significantly during the measurement. Alternatively, the electrolyte can be pumped with electrolyte pumps.

The measuring cells may only be operated in a laboratory environment.

## Safety instructions

The test cell ElyFlow is only intended for the above-mentioned purpose.

The manufacturer is not liable for damage resulting from improper use.

Observe the following instructions, otherwise the test cell may be damaged, or measurement results may be falsified.

When working with chemicals, all relevant safety regulations of the manufacturer and the laboratory must be observed.

The cell heater and / or the cell current you apply can heat the test cell over 80°C.

The test cell must be placed on a base that is insensitive to heat.

Do not place heat sensitive things on the test cell.

Do not place inflammable things near the test cell.

Pay attention to a correct lining of your measurement and heating cable.

Select cable insulations with a reasonable temperature specification. Else there is a risk of an electrical shortcut and fire hazard.

Electrolyte vapors are constantly emitted from the test cell according to the vapour pressure curve.

The test cell should be operated on a laboratory tray that can safely hold the entire filled electrolyte volume in case of leakage.

Protect electrical / electronic equipment by a large distance, partitions, or operation under a fume cupboard.

Wear safety glasses.

Wear protective gloves.


The operator of the test cell must be equipped with adequate laboratory protective equipment according to the safety data sheets and specifications for your electrolyte.

## Use

### Device overview

#### Layout

## Electrochemical test cell ElyFlow

In contrast to the usual test cells, Gaskatel test cell ElyFlow is made of PTFE, because glass is a poor material. Firstly, it breaks quickly. Secondly, it is chemically stable only to a limited extent because silicates dissolve. This is a particular problem with aggressive solutions such as alkalis.

But also for light sensitive electrolytes glass is an insufficient material.

In very poorly conductive electrolytes, the conductive silicate layers cause an uneven field line distribution.

Like the Gaskatel test cells FlexCell, ElyFlow is also characterised by high mechanical stability and high chemical resistance. The plastics have no influence on the field line distribution.


In ElyFlow, as in the FlexCell measuring cells, the electrolyte chamber has been designed in a tubular shape.

ElyFlow is also designed so that you can pump your electrolyte from a storage vessel through the measuring cell with the help of a diaphragm pump.

This allows longer measurement times. The flow rate should be less than 0.1 l/min.

In addition, the electrolyte in the storage vessel can be continuously monitored during the measurement: Temperature, conductivity, pH value, analysis of reaction products.

Use a heatable storage vessel when measuring at higher temperatures.

### Gaskatel test cells in comparison

The test cell ElyFlow with 10 cm<sup>2</sup> has got a larger working surface than FlexCell test cells. The distance between working and counter electrode is in ElyFlow 20 mm and therefore much smaller than in FlexCell. This allows measurements in less conductive electrolytes.

	ElyFlow	FlexCell
Parallel three electrode set-up	✓	✓
Integrated Haber-Luggin-capillaries	✓	✓
Controlled remove of reaction gases	✓	✓
Distance between working and counter electrode	20 mm	70 mm
Working surface	10 cm <sup>2</sup>	3 cm <sup>2</sup>
Sample size	65 mm x 70 mm	30 mm x 50 mm
Electrolyte volume	25 ml	30 ml
Counter electrode	Nickel; MMO	Platinum/Iridium
Counter electrode replaceable	✓	(✓)
Can be temperature controlled	✓	✓
Gas supply for gas diffusion electrodes	✓	✓
Measurements on gas diffusion electrodes	✓	✓
Measurements on plates	✓	✓
Electrolyte circuit possible	✓	✗
Maximum recommended measurements duration	Unlimited with electrolyte circuit	24 h
Measurements on membranes with additional analyte compartment	✓	(✓)
Determination of hydrogen permeation with additional analyte compartment	✓	✗

## Test cell ElyFlow in detail

### Working electrode


The working electrode is placed between working electrode holder and analyte compartment. The seal is made via the O-rings inserted in both parts. Alternatively suitable silicone seals can also be used.

The sample must be cut to the dimensions 65 mm x 70 mm. Then the sample fits exactly between the fixing pins that prevent the sample from slipping. If you want to use electrodes of other geometries, you can simply remove these pins.

If the working electrode is a gas diffusion electrode with a non-conductive gas diffusion layer such as porous PTFE, then the PTFE foil must be removed approx. 1 cm from the edge so that the sample can be contacted there.

### Gas supply


For measurements on gas diffusion electrodes, these must be supplied with gas. The working electrode holder has two connections: the gas inlet at the top and the gas outlet at the side. If liquid gets through the electrodes, it can be removed with the gas flow.


### Contact working electrode

Two gold-plated contact screws (4 mm multilam plugs) press onto the working electrode from the working electrode holder's backside. If the contact is not sufficient, the multilam plugs can be tightened with a socket spanner.

Make sure that these plugs are clean. Replace these plugs if the gold-plated layer no longer looks sufficiently good. This also applies to all test leads.

Generally, the working electrode holder is made of PTFE. Optional Gaskatel offers a transparent working electrode holder made of PMMA (plexiglass). This one allows a better view on the working electrode and is recommended for measurements on gas diffusion electrodes.

### Analyte compartment

Only test cells with a tubular arrangement with equally sized measuring and counter electrodes guarantee a parallel field line distribution.


Parallel field lines in turn guarantee parallel equipotential surfaces and thus a reproducible potential measurement, as the reference electrode measures on these potential surfaces.

Gaskatel test cells have such parallel field lines, as the working and counter electrodes have the same size and are arranged parallel to each other due to the tubular electrolyte space.


The analyte compartment is tightened with an O-ring towards the counter electrode.

Another O-ring seals towards the sample. Fixing pins ensure that the sample is held in place during assembly.


The fixing pins can also be inserted in the analyte compartment instead of the working electrode holder. They provide a hold for the sample during assembly.

Furthermore, there is an electrolyte inlet and outlet to allow electrolyte to flow through the cell.

The electrolyte inlet is divided into a total of six small nozzles. Three point to the counter electrode, the other three to the working electrode. The electrolyte outlet is located at the top of the cell and is provided with a cavity (dome) that ensures the removal of gas bubbles originating from the electro-chemical reaction inside the cell. Where glass bubbles adhere, there can be no electrolyte. The reaction surface is shadowed. Local elements can develop. Three Haber-Luggin capillaries are integrated at half height in the immediate vicinity of the working electrode.

### **Reference electrode**

Gaskatel test cells are optimized for the use of the reversible hydrogen reference electrodes (RHE) Mini-HydroFlex.

In contrast to the usually used calomel or silver electrodes, these hydrogen reference electrodes do not introduce any foreign electrolyte into the measuring solution. This is very important for corrosion measurements, for example, where the smallest amounts of chloride can drastically influence and falsify the results.

Most electrochemical reactions are pH dependent. This is compensated using Mini-HydroFlex, because Mini-HydroFlex itself is also pH-dependent.

### **Position of the reference electrode – Haber-Luggin-capillary**

For precise potential measurements it is necessary that the reference electrode is positioned close to the working electrode.

To reduce the voltage drop across the electrolyte, the so-called IR drop, Haber-Luggin capillaries are used. However, for measurement results to be comparable between cells, the Haber-Luggin capillary must meet very precise and reproducible requirements.

The field lines in the electrolyte must not be disturbed (small diameter of the Haber-Luggin-capillary). In turn, a large diameter of the Haber-Luggin-capillary is necessary for low electrolyte resistance. Therefore, the Haber-Luggin-capillaries in Gaskatel test cells are located outside the electrolyte space. This allows a slightly larger diameter without disturbing the field lines.

The capillary must end at a defined, fixed distance from the working electrode. From a quality assurance point of view, it is also important that all these factors are kept constant for different half-cells.

This is where the precision of CNC manufacturing in plastic helps. The Gaskatel measuring cells have very precise Haber-Luggin capillaries that end at the working electrode. This makes the measurements very precise and comparable.


The capillaries must not be blocked by gas bubbles, as this leads to contact losses during the potential measurement. The resistance then increases dramatically and the potentiostat starts to oscillate.

Therefore, the Haber-Luggin capillaries in the FlexCell measuring cells are filled with a solid electrolyte.

### Counter electrode

ElyFlow is supplied with a counter electrode made of nickel and a counter electrode with MMO coating on titanium. However, you can use other counter electrodes with the dimensions 64 mm x 64 mm x 2 mm.

The contact is made on the back via two gold-plated 4 mm multilam plugs.

### Temperature

Many processes take place at higher temperatures. Therefore, it is often necessary to heat the test cells.

Gaskatel ElyFlow test cell can be heated by a heating element, integrated behind the counter electrode.

Please use our temperature control box 83900 for temperature control or the control box TP 83700, because only a well-regulated temperature ensures comparable measurement results.

### Crevice corrosion

If the sample is not sealed properly, an electrochemical reaction may take place underneath the sealing or even outside of the active area. A proper sealing is essential. Oxygen may only enter in a controlled manner.

The measuring cells are therefore supplied with the appropriate O-rings.


Silicone sealings are also suitable and available cut to size.

### **Reaction products of the counter electrode – additional analyte compartment**

Often very aggressive ions are produced at the counter electrode. Depending on the electrolyte, peroxides, perchlorates or persulphates etc. can be produced. Even in the smallest quantities may initiate corrosion on the working electrode.

An additional analyte compartment (83420) provides a remedy here and can be used to install an additional membrane between the counter electrode and working electrode holder. In this way, the working and counter electrodes are separated from each other.

### **Control box**

Gaskatel has developed its own control box (control box TP 83700). It can be used to control the temperature of the cell.

An integrated galvanostat and voltmeter allow to perform simple electrochemical measurements.

Additionally the control box has two control units for electrolyte pumps.

The complete cell control is supplied with a USB connection to your computer.

A corresponding control software is available for download.

### **Electrolyte pumps**

Pumping hot and often corrosive electrolyte in lab size is not an easy way. There are some commercial pumps available you may use: impeller pumps, hose pumps and membrane pumps. When using hose pumps, you should use EDPM or FKM tubes. In the case of membrane pumps, the preferred materials are PP or PTFE for the housing and EDPM or FKM as a membrane material.

Gaskatel is offering a membrane pump out of PP housing and EPDM membrane with an internal regulated power supply (83721).

## Operation

### 1. Assembling the test cell

First, the working electrode may have to be cut to size:


Width 6,5 cm, length 7 cm

The thickness of the working electrode is almost arbitrary.

The active area in the electrolyte is 10 cm<sup>2</sup>.

The contact screws in the working electrode holder must be turned out a little so that the electrode cannot touch these contacts at first.

The cut working electrode is placed between the two silicone sealings or EPDM O-rings. The fixing pins ensure that the electrode cannot slip.

When screwing the cell components together, make sure that the screws are tightened diagonally to each other one after the other. There must be no gap between the components. Tighten the screws only hand-tight, especially when using the flat gaskets. Leakage is usually caused by unevenly tightened screws and not by too weakly tightened ones.


If this is the case, you can now screw the contact screws back in with the socket spanner and carefully tighten them so that you contact the electrode.

Caution is required here, especially with very thin samples.


## 2. Use of the analyte compartment with membrane

In addition, an intermediate plate (analyte compartment) can be inserted, which one is optionally available and supplied with the necessary longer screws.

A membrane is inserted between the analyte compartment and the counter electrolyte compartment. This must be cut round and 50 mm in diameter. This membrane separates anolyte and catholyte. The working electrode is thus protected from the reaction products of the counter electrode. When selecting the membrane, pay attention to the electrolyte resistance and the ion-exchanging properties of the membrane.


You can also measure membrane resistances with this arrangement. Two reference electrodes are then needed. One is inserted into the analyte compartment and the second into the bore of the electrolyte compartment as usual.


### 3. Filling of the test cell

ElyFlow is designed to operate with electrolyte pumps. Therefore, all openings must be closed.


Before filling the cell, the Mini-HydroFlex must be screwed in into the reservoir marked with RE1. If a temperature sensor is to be used, this must also be screwed in, now. For each analyte compartment an individual pump is necessary. The connection is made via 6/4 mm hoses, fitting to the supplied straight connectors.

If the hoses are connected safely at cell and pumps, if the suction and the outlet hoses are hanging in the corresponding vessels, the pumps can be started, with low rate to flood the single compartments carefully. In case you wish to control the level in the reservoirs for Mini-HydroFlex, switch-off the pumps. screw-out the Mini-HydroFlex a little, so you can look inside the reservoir.

Without the pumps, the straight connector for the inlet must be replaced by a screw plug. After removing the outlet straight connector in the center opening, the cell can be filled with approx. 25 ml of electrolyte. Please consider: the reservoirs for the Mini-HydroFlex are filling by themselves, but with a delay in time. They could be filled half to avoid overflowing when setting the reference electrode.

If the hole does not fill with electrolyte, air bubbles in the Haber-Luggin capillary block the electrolyte flow. Do not insert the reference electrode before the reservoirs are half filled.

#### 4. Cleaning of the test cell

Usually, it is sufficient to rinse the cell thoroughly with water after disassembly. Deposits can be removed mechanically with a suitable laboratory brush. Stubborn metal deposits such as silver or copper can be removed with nitric acid. In the analyte compartments, the O-rings and the grooves in which they are located should be thoroughly rinsed and cleaned after the measurements.

#### 5. Use of the ElyFlow-System

##### Characterisation of electrodes

With the ElyFlow system you can perform electrochemical measurements on samples such as metals or gas diffusion electrodes. Since it is designed to pump electrolytes through, long-term measurements are possible with ElyFlow.


##### Measure membrane resistance

ElyFlow is suitable for measuring membrane resistance. Two analyte compartments and two reference electrodes are required for this.

An electrode (e.g. metal sheet) must be installed between the working electrode holder and analyte compartment 1.


The membrane is installed between the two analyte chambers. It must be cut round and 50 mm in diameter.

The nickel sheet in the counter electrode holder acts as a second electrode.

Now it is important that the reference electrodes are inserted correctly. They must be inserted in the reservoirs where the Haber-Luggin capillaries face the membrane. One reference electrode must be inserted in the reservoir with the sticker RE2 in analyte compartment 1 as the Haber-Luggin capillary is aligned there towards the membrane.

The second reference electrode must be inserted in the reservoir with the sticker RE1 in analyte compartment 2.

### Galvanic coatings

Electroplating of your samples can also be carried out with the ElyFlow system, as you can assemble the system to suit your task by combining the individual modules.

### Measuring hydrogen permeation

The ElyFlow-system is excellent for measuring hydrogen permeation.

Together with **IChemAnalytics GmbH**, Gaskatel offers a complete measurement system, consisting of the ElyFlow test cell from Gaskatel and the Automated Measurement and Control Box (AMB) from **IChemAnalytics GmbH**.

When examining steel foils, it is also possible without Pd coating on the oxidation side of the sample.

With this system, you can perform a wide range of electrochemical measurement and control tasks:


- In situ monitoring of hydrogen evolution during electrochemical coating, including monitoring of the potential at the phase boundary.
- In situ monitoring of hydrogen evolution in pre-treatment processes, such as pickling processes, in conjunction with monitoring of the electrochemical potential at the phase boundary.
- Monitoring of production process parameters with regards to hydrogen embrittlement.
- Determination of H-diffusion coefficients in metallic materials in dependence of the morphological material properties.

Further information you will find at [www.gaskatel.com](http://www.gaskatel.com).

### Set-up for a fuel cell

The individual modules of the ElyFlow system can also be combined to build a fuel cell. To do this, two working electrode holders and one analyte compartments must be combined. Working electrode holders made of PMMA allow a view of the gas diffusion electrodes. For example, it can be observed whether liquid drops form on the gas side of the electrodes because the electrolyte penetrates (so-called "flooding" of the electrode) or whether water cannot be removed.

Temperature control via the cell components is not possible in this set-up. The electrolyte can then be tempered in a storage vessel and must then be pumped through the cell. The temperature should not exceed 80°C because of the PMMA components.

## 6. Troubleshooting

Error	Possible causes	Check
Cell temperature does not exceed	Poor contact	Check the connecting leads of the heating elements
	PTC heating element in silicone is defect	Measure the resistance of cold the heating elements: 15 Ohm
Potentiostat is going in overpotential	Poor contacts	Check the connecting leads of the potentiostat Check the multilam plugs of the cell
	Counter electrode	Remove the counter electrode and check its condition
False potential	Gas bubbles in front of the reference electrode	Move the reference electrode a little up and down
	Gas bubbles in Haber-Luggin capillary	Suck air out of the reference electrode hole with a pipette ball
Noisy signal	Impedance of the reference electrode or Haber-Luggin capillary is too high	Adjust the current range of the potentiostat
Electrolyte leakage	Plastic thumb screws not tightened	Tighten the thumb screws
	O-ring or silicone sealing is not placed properly	Check the position of the contact screws and fixing pins.

## Operating conditions

Nominal voltage working- / counter electrode:	< 50 VAC/< 75 VDC
Nominal current working- / counter electrode:	max. 19 A
pH-range:	pH -2 bis pH 16
Temperature:	10°C bis 100°C
Pressure:	max. 200 mbar

## Technical data

Dimensions: without analyte compartment, protruding screw connections, screws	124 mm x 80 mm x 95 mm
Electrolyte volume:	25 ml
Material counter electrode holder:	Polytetrafluorethylene (PTFE)
Material working electrode holder:	Polytetrafluorethylene (PTFE)
Material analyte compartment	Polytetrafluorethylene (PTFE)
Active sample area:	10 cm <sup>2</sup>
Working electrode sample size:	6,5 x 7 cm
Contacts working electrode:	2 multilam plugs 4 mm, gold-plated
Sealings working electrode	O-Rings EPDM
Screws fitting gas supply:	ECTFE G 1/8" for 6/4 mm hose
Heating element	70W, 30V
Counter electrode:	Nickel; Mixed metal oxide (MMO)
Contacts counter electrode	2 multilam plugs 4 mm, gold-plated
Fixing screws	Polyamide PA 6.6

## Disposal

The test cell must not be disposed of with household waste. Please dispose of the test cell in accordance with local regulations in a separate collection for electrical and electronic equipment.

If sent to us free of charge (sufficiently stamped package), we will dispose of the used test cell from our product range free of charge.

## Exclusion of liability

The product is subject to continuous further development. For this reason, it is possible that parts of the instructions, technical data or images in these instructions may differ slightly from the product in front of you. The information in these instructions is only intended to clarify how the product is to be handled.

If you have any questions, please do not hesitate to contact us.

Legal claims based on these instructions cannot be asserted.