
Mini-HydroFlex

Wasserstoffreferenzelektrode
Hydrogen Reference Electrode

Artikelnummer / Item number 81020 (81021)

User Manual

Version 2.3

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Lilienthalstrasse 146

D-34123 Kassel

+49 (0) 561 59190

info@gaskatel.de

www.gaskatel.com

Languages

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Operating instructions

Table of contents

| | |
|---------------------------------|----|
| Overview | 2 |
| Delivery scope | 2 |
| Declaration of conformity | 2 |
| General information | 2 |
| About these instructions | 2 |
| Intended use | 2 |
| Safety instructions | 3 |
| Use | 3 |
| Device overview | 3 |
| Operation | 4 |
| Operating conditions | 10 |
| Technical data | 10 |
| Disposal | 11 |
| Exclusion of liability | 11 |

Overview

Delivery scope

- 1 x Mini-Hydrogen reference electrode inclusive
hydrogen cell (80521) and activation sticker (80077)
- 1 x Manual

Declaration of conformity

 This product complies with all EU directives applicable to this product.

General information

About these instructions

Read these instructions carefully. To ensure a long service life and reliable use, all instructions mentioned in this manual must be observed and followed. Please keep these instructions at hand.

You can access the instructions and further information at www.gaskatel.com.

Intended use

Please consider the Good Laboratory Practice (GLP).

Mini-HydroFlex is a reference electrode, exclusively intended for potential measurement during electrochemical measurements.

Safety instructions

The electrode is only intended for the above-mentioned purpose. The manufacturer is not liable for damage resulting from improper use.


Respect the following instructions, otherwise the electrode may be damaged, or the measurement results may be falsified.

When working with chemicals, all relevant safety regulations of the manufacturer and the laboratory must be respected.

Use

Device overview

Construction and dimensions

Properties of the hydrogen reference electrode Mini-HydroFlex

The hydrogen reference electrode is made entirely from plastic. Therefore, it is robust and insensitive to breakage.

It can be used in **highly concentrated alkaline solutions** as well as in **concentrated acids**, even in **hydrofluoric acid or fluoride-containing media**.


Of course, Mini-HydroFlex is **free of mercury and silver**.

Operation

1. Commissioning

Unscrew the hydrogen cell holder and remove the hydrogen cell.

Take the replacement hydrogen cell and remove the protective foil. There are two openings visible. If any adhesive residue remains, remove it by rubbing with a cloth.

Insert the hydrogen cell with these two openings down into the recess.

Turn the hydrogen cell holder back in.

Mark month and year on the activation sticker. Attach the sticker to the electrode head.

Remove the protective cap from the electrode carefully.

Place Mini-HydroFlex in distilled water for 24 hours.

Please do not remove the labelling with the serial number and the activation sticker with the activation date.

2. Storage of Hydrogen Reference Electrode Mini-HydroFlex

After the measurements, rinse the electrode thoroughly with water. Place the hydrogen electrode in a liquid (such as measuring solution, 1 mol/l hydrochloric acid, 1 mol/l sodium hydroxide solution, water) to keep it functional.

As soon as the runtime of 12 months is reached, the hydrogen cell must be replaced immediately (see 4. Maintenance). Otherwise, a negative pressure will build up in the capillary so that liquid can get into the capillary.

3. Checking Mini-HydroFlex

If you want to check Mini-HydroFlex, we recommend the standard calomel electrode or, if necessary, the silver silver chloride electrode.

We recommend the following solutions as measuring solution:

| | | |
|---|------|--------|
| Hydrochlorid acid | pH 0 | 244 mV |
| Hydrochloric acid/citric acid | pH 2 | 362 mV |
| Citric acid/sodium chloride/sodium hydroxide solution | pH 4 | 480 mV |

The potentials given refer to the standard calomel electrode at 25 °C.

4. Maintenance


Replacement of the Hydrogen cell

Unscrew the hydrogen cell holder and remove the hydrogen cell.

In the recess in the electrode head, you see a small O-ring. Please grease this, e. g. with silicone grease.

Take the replacement hydrogen cell and remove the protective foil. There are two openings visible. If any adhesive residue remains, remove it by rubbing with a cloth.

Insert the hydrogen cell with these two openings down into the recess.

Turn the hydrogen cell holder back in.

Mark month and year on the activation sticker. Attach the sticker to the electrode head.

Place Mini-HydroFlex in distilled water for 24 hours.

5. Troubleshooting

Potential is not displayed correctly

Cause 1: poor contacting

Check the test leads on optical damages such as corrosion, cracks, and sessile plugs. Replace the test leads.

Check the electrode with an external reference electrode, e.g. a calomel or silver silver chloride electrode (see 0. Checking Mini-HydroFlex).

Cause 2: no or too little hydrogen production

Did you place the new Mini-HydroFlex in water for 24 hours?

If not, wait 24 hours before you measure.

Is the duration of the hydrogen cell exceeded?

If yes, please replace the hydrogen cell (see 4. Maintenance).

Cause 3: slow ion exchange, e.g. when you switch from concentrated solutions to weak concentrated solutions

Wait the setting time, sometimes the balancing of concentrations lasts longer than expected. If necessary, check the potential in another electrolyte such as 1 mol/l hydrochloric acid (see 0. Checking Mini-HydroFlex).

Cause 4: air/oxygen gets to the electrode

Avoid that gases such as air or oxygen reach the hydrogen electrode's shaft. If this happens the hydrogen is displaced or reacts off. No hydrogen potential can be set.

Fluctuating, noisy, or oscillating potentials

Cause 1: hydrogen bubbles of the hydrogen electrode itself

Bubbles are continuously leaking out of the hydrogen electrode. Some bubbles are small, others are bigger. Usually, they do not disturb your measurements.

Is a big bubble forming at the bottom of the hydrogen electrode which sticks to the vessel wall?

If possible, place the electrode further away from the vessel wall or hang the hydrogen electrode transversely into the measuring vessel.

Cause 2: induced gases

Change the position of your gas inlet.

Induced gas bubbles which are piped by close to the hydrogen electrode can disturb the measurements and lead to fluctuating potentials.

Cause 3: potentiostat or measuring device

Check your measuring device.

If you measure in low conductive electrolytes, potentiostats and measuring devices quickly reach their limits.

Deposits (e.g. red or brown) on the electrode

The hydrogen electrode is contaminated, and a mixed potential is set

Check the runtime of the hydrogen cell. If necessary, the hydrogen cell must be replaced as described in 4. Maintenance.

Some deposits can be wiped off with a cloth.

With 1000 grit sandpaper, you can carefully sand the bottom of the electrode.

More stubborn impurities dissolve in 1 mol/l nitric acid. The electrode must then be rinsed thoroughly in distilled water.

If that is not successful, the electrode is so badly contaminated that we suggest a new purchase.

6. Further information

For more information, please visit www.gaskatel.com.

Operating conditions

| | |
|--|-------------------|
| pH range: | pH -2 to pH 16 |
| Temperature: | -20 °C to 120 °C |
| HF/Fluoride: | Up to 500 000 ppm |
| Caustic potash, caustic soda solution: | Up to 50 wt% |

Technical data

| | |
|------------------------|-------------------------------|
| Total length: | 87 ± 5 mm |
| Shaft length: | 53,5 ± 0,5 mm |
| Shaft diameter: | 1.6 mm |
| Electrode shaft: | Polyether ether ketone (PEEK) |
| Electrode head: | Polytetrafluorethylene (PTFE) |
| Thread: | G 1/8" |
| Measuring electrode: | Platinum/Palladium/PTFE |
| Hydrogen cell holder: | Polyvinyl chloride (PVC) |
| Runtime hydrogen cell: | 12 months |

Disposal

The electrode must not be disposed of with household waste. Please dispose the electrode in accordance with local regulations in a separate collection for electrical and electronic equipment.

If sent to us free of charge (sufficiently stamped package), we will dispose of the used electrodes from our product range free of charge.

Exclusion of liability

The product is subject to continuous further development. For this reason, it is possible that parts of the instructions, technical data or images in these instructions may differ slightly from the product in front of you. The information in these instructions is only intended to clarify how the product is to be handled.

If you have any questions, please do not hesitate to contact us.

Legal claims based on these instructions cannot be asserted.