
# pHydrunio

Wasserstoff-pH-Elektrode  
Hydrogen pH Electrode

Artikelnummer / Item number 84100

User Manual

Version 2.0

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Lilienthalstrasse 146

D-34121 Kassel

+49 (0) 561 59190

info@gaskatel.de

[www.gaskatel.com](http://www.gaskatel.com)

## Languages

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# Operating Instructions

## Table of contents

Overview .....	2
Delivery scope .....	2
Declaration of conformity .....	2
General.....	2
About these instructions .....	2
Intended use .....	2
Safety instructions.....	2
Use .....	3
Device overview.....	3
Operation .....	3
Operating conditions .....	11
Technical data.....	11
Disposal.....	12
Exclusion of liability .....	12

## Overview

### Delivery scope

- 1 pHDrunio Hydrogen pH electrode inclusive hydrogen cells (80521)
- 1 Connection lead (89015)
- 100 ml Inner buffer (84210)
- 1 Syringe 10 ml with cannula (80341)
- 1 Pressure control valve (80189)
- 1 Manual

### Declaration of conformity

 This product complies with all EU directives applicable to this product.  
pHDrunio has a food approval.

## General

### About these instructions

Read these instructions carefully. To ensure a long service life and reliable use, all instructions mentioned in this manual must be observed and followed. Please keep these instructions at hand.

You can access the instructions and further information at [www.gaskatel.com](http://www.gaskatel.com).

### Intended use

Please consider the Good Laboratory Practice (GLP).

The hydrogen pH electrode is intended exclusively for measuring pH values in a pH range from -2 to 16.

## Safety instructions

The electrode is only intended for the above-mentioned purpose. The manufacturer is not liable for damage resulting from improper use.


Respect the following instructions, otherwise the electrode may be damaged, or the measurement results may be falsified.

When working with chemicals, all relevant safety regulations of the manufacturer and the laboratory must be respected.

## Use

### Device overview


#### Construction and dimensions

## Principle of the hydrogen pH electrode pHydrunio

pHydrunio is the world's first commercially available platinum hydrogen pH electrode (patent DE 10 2011 113941). It combines a reversible platinum hydrogen electrode as a measuring electrode with a platinum hydrogen reference electrode in an internal buffer as a reference electrode.

The platinum hydrogen electrode is the only electrode whose potential is directly determined by the hydrogen ion activity. pH value and hydrogen potential are directly linked. The potential, which depends on the pH value, can be calculated with the help of Nernst equation.

$$E = E^0 + \frac{RT}{2F} \ln \left[ \frac{a^2(H^+)}{p(H_2)} \right]$$

$$E = \frac{RT}{F} \ln[a(H^+)]$$

$$E = 2.303 \frac{RT}{F} \log[a(H^+)]$$

$$E = -2.303 \frac{RT}{F} pH$$

$$E = 0.058 \text{ V @ } 20^\circ\text{C} \quad E = 0.059 \text{ V @ } 25^\circ\text{C} \quad E = 0.060 \text{ V @ } 30^\circ\text{C}$$

Measuring voltage U and pH value are related to each other as follows:

$$U = -2.303 \frac{RT}{F} (pH_{inside} - pH)$$

The pH value can be calculated from the measuring voltage by conversion.

$$pH = 7 - \frac{U}{0.059 \text{ V}} \quad @25^\circ\text{C}$$

*E* = Hydrogen potential in Volt

*E*<sup>0</sup> = Standard Hydrogen potential in Volt = 0,0 V per Definition

*R* = 8.314 K mol<sup>-1</sup>K<sup>-1</sup>

*T* = Temperature in K

*F* = 96485 C mol<sup>-1</sup>

*p*(H<sub>2</sub>) = Saturated vapour pressure = 1000 hPa

*a* = Activity in mol/l

*U* = Measuring voltage pH-Electrode in Volt

*pH*<sub>inside</sub> = 7

## Properties of the hydrogen pH electrode pHydrunio

The pH electrode pHydrunio is made entirely of plastic. Therefore, it's robust and insensitive to breakage.

It can be used in **highly concentrated alkaline solutions** as well as in **concentrated acids**, even in **hydrofluoric acid or fluoride-containing** media.

Since it is **glass-free**, it shows **neither an acid nor an alkali error**.

In contrast to conventional pH glass electrodes, it shows **no phosphate error**.

Of course, pHydrunio is **free of mercury and silver**.

## Operation

### 1. Commissioning

The pHydrunio platinum hydrogen pH electrode is supplied pre-assembled, i.e. the hydrogen cells are already installed and activated. There is a sticker on the electrode head indicating the month of activation. This must not be removed.

pHydrunio is already filled with internal electrolyte. Replace the red transport plug in the electrode head with the enclosed red pressure relief valve, which is also a sealing plug.

Carefully remove the transport cap from the bottom of the electrode shaft.

Place the hydrogen pH electrode in distilled water for 24 hours.

## 2. Storage of the hydrogen pH electrode

After the measurements, rinse the hydrogen pH electrode thoroughly with water. Please place the pH electrode in a vessel with distilled water. Under no circumstances it should be stored in the transport cap. The hydrogen would squeeze slowly the small volume of liquid out of the cap. The hydrogen pH electrode should not be stored dry in air as long as the hydrogen cells are generating hydrogen. A mixed potential forms which can lead to incorrect results.


As soon as the running time of 12 months is reached, the hydrogen cells must be replaced immediately (see 5. Maintenance). Otherwise, a negative pressure will build up in the capillaries so that liquid can get into the capillaries.

Please do not remove the labelling with the serial number and the sticker with the activation date.

## 3. Calibration

pHydrunio can be used in most common buffers. Note the slightly longer response times of our pH electrode.

#### 4. Use with pH meters and measuring instruments

pHydrunio can be connected to any pH meter with analogue input (BNC or open wire ends).

In contrast to the usual glass electrodes, our hydrogen pH electrode has a low resistance.

#### 5. Maintenance

##### Replacement of hydrogen cells


First, please unscrew the red cap on the electrode head.

There are two black caps on the sides of the electrode head. Please also unscrew these. Take out the hydrogen cells underneath.

You will see a small O-ring in the recess. Please grease this, e.g. with silicone grease.

Take a new hydrogen cell and remove the protective foil. Two openings become visible. If any adhesive residue remains, remove it by rubbing with a cloth.

Place one hydrogen cell in the recess with these openings facing down.

Screw the black caps back in.

Turn the red pressure relief valve back into the electrode head.

Place the hydrogen pH electrode in distilled water for 24 hours.

Mark month and year on the label. Attach the sticker to the electrode head.

### Changing the inner buffer

The internal buffer (reference electrolyte) of pHDrunio can be replaced. To do this, unscrew the red pressure relief valve on the electrode head. You can pull out the internal buffer with the enclosed syringe and cannula. Rinse the electrode with distilled water. Fill in new inner buffer directly from the vial.

Carefully screw the pHDrunio back on.

Place pHDrunio in distilled water for 24 hours.

## 6. Troubleshooting

### **pH value is not displayed correctly**

#### Cause 1: buffer solution

Please measure in fresh buffer solution.

If buffer solutions are exposed to air for a longer time, the pH value can shift.

This occurs mainly with alkaline buffer solutions.

#### Cause 2: slow ion exchange, e.g. when changing from concentrated solutions to weakly concentrated solutions

Await the setting time, sometimes the balancing of concentrations lasts longer than expected. If necessary, check the potential in another electrolyte, e.g. in calibration buffers.

#### Cause 3: internal buffer contaminated

It can happen that measuring solution enters the internal buffer and contaminates the internal buffer.

In that case, the internal buffer has to be replaced (see 5. Maintenance).

#### Cause 4: bad contacting

Check the measuring cables on optical damages such as corrosion, cracks and sessile plugs. Replace the cable. Check both electrodes with an external reference electrode, e. g. HydroFlex, calomel or silver silver chloride electrode.

#### Cause 5: air/oxygen gets to the hydrogen electrode

Avoid that gases such as air or oxygen reach the hydrogen electrode 's shaft.

If this happens, the hydrogen is displaced or reacts off. No hydrogen potential can be set.

#### Cause 6: no or too low hydrogen production

Is the runtime of the hydrogen cells exceeded?

If yes, please replace the hydrogen cells (see 5. Maintenance).

## **Fluctuating, noisy or oscillating pH values (potentials)**

### Cause 1: hydrogen bubbles of the hydrogen electrode itself

Bubbles are continuously leaking out of the platinum hydrogen electrode. Some bubbles are small, others are bigger. Usually, they do not disturb your measurements.

Is a big bubble forming at the bottom of the hydrogen electrode which sticks to the vessel wall? If Possible, place the electrode further away from the vessel wall or hang the hydrogen electrode transversely into the measuring vessel.

### Cause 2: induced gases

Change the position of your gas inlet. Induced gas bubbles which are piped by close to the hydrogen electrode can disturb the measurements and lead to fluctuating potentials.

## **7. Further information**

For more information, please visit [www.gaskatel.com](http://www.gaskatel.com).

## Operating conditions

pH range:	pH -2 to pH 16
Temperature:	-20 °C to 120 °C
Pressure:	Low overpressure possible
HF/Fluoride:	Up to 500 000 ppm
Caustic potash, caustic soda solution:	Up to 50 wt%

## Technical data

Total length:	185 mm
Shaft length:	120 mm
Shaft diameter:	12.2 mm
Diaphragm:	PTFE/ceramic
Inner buffer:	KCl/Buffer based on phosphate
Electrode shaft:	Polypropylene
Electrode head:	Polypropylene
Thread:	PG 13.5
Measuring electrode, reference electrode:	Platinum, Palladium, PTFE
Runtime hydrogen cells:	12 Month

## Disposal

The electrode must not be disposed of with household waste. Please dispose of the electrode in accordance with local regulations in a separate collection for electrical and electronic equipment.

If sent to us free of charge (sufficiently stamped package), we will dispose of the used electrodes from our product range free of charge.

## Exclusion of liability

The product is subject to continuous further development. For this reason, it is possible that parts of the instructions, technical data or images in these instructions may differ slightly from the product in front of you. The information in these instructions is only intended to clarify how the product is to be handled.

If you have any questions, please do not hesitate to contact us.

Legal claims based on these instructions cannot be asserted.