
Mini-HydroFlex

Wasserstoffreferenzelektrode
Hydrogen Reference Electrode

Artikelnummer / Item number 81020 / 81021

User Manual

Version 2.1

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Lilienthalstrasse 146

D-34123 Kassel

+49 (0) 561 59190

info@gaskatel.de

www.gaskatel.com

Languages

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Operating Instructions

Table of contents

Overview	2
Delivery scope	2
Declaration of conformity	2
General.....	2
About these instructions	2
Intended use	2
Safety instructions.....	2
Use	3
Device overview	3
Operation	3
Operating conditions	9
Technical data.....	9
Disposal.....	10
Exclusion of liability	10

Overview

Delivery scope

- 1 Mini-Hydrogen reference electrode inclusive hydrogen cell (80521)
- 1 Manual

Declaration of conformity

 This product complies with all EU directives applicable to this product.

General

About these instructions

Read these instructions carefully. To ensure a long service life and reliable use, all instructions mentioned in this manual must be observed and followed. Please keep these instructions handy.

You can access the instructions and further information at www.gaskatel.com.

Intended use

Please consider the Good Laboratory Practice (GLP).

The Mini-HydroFlex is a reference electrode, exclusively intended for potential measurement during electrochemical measurements.

Safety instructions

The electrode is only intended for the above-mentioned purpose. The manufacturer is not liable for damage resulting from improper use.


Respect the following instructions, otherwise the electrode may be damaged or the measurement results may be falsified.

When working with chemicals, all relevant safety regulations of the manufacturer and the laboratory must be respected.

Use

Device overview

Construction and dimensions

Properties of the Hydrogen Reference Electrode Mini-HydroFlex

The hydrogen reference electrode is made entirely from plastic. Therefore, it's robust and insensitive to breakage.

It can be used in **highly concentrated alkaline solutions** as well as in **concentrated acids**, even in **hydrofluoric acid or fluoride-containing media**.


Of course, Mini-HydroFlex is **free of mercury and silver**.

Operation

1. Commissioning

Unscrew the screw cap and remove the hydrogen cell.

Take the replacement hydrogen cell and remove the protective foil. There are two openings visible. If any adhesive residue remains, remove it by rubbing with a cloth.

Insert the hydrogen cell with these two openings down into the recess.

Turn the screw cap back in.

Mark month and year on the label. Attach the sticker to the electrode head.

Remove the protective cap from the electrode carefully.

Place Mini-HydroFlex in water for 24 hours.

Please do not remove the labelling with the serial number and the sticker with the activation date.

2. Storage of Hydrogen Reference Electrode Mini-HydroFlex

After the measurements, rinse the electrode thoroughly with water. Place the hydrogen electrode in a liquid (such as measuring solution, 1 mol / l hydrochloric acid, 1 mol / l sodium hydroxide solution, water) to keep it functional.

As soon as the running time of 12 months is reached, the hydrogen cell must be replaced immediately (see 4. Maintenance). Otherwise, a negative pressure will build up in the capillary so that liquid can get into the capillary.

3. Checking the Mini-HydroFlex

If you want to check HydroFlex, we recommend the standard calomel electrode or, if necessary, the silver silver chloride electrode.

We recommend the following solutions as measuring solution:

Hydrochlorid acid	pH 0	244 mV
Hydrochloric acid/citric acid	pH 2	362 mV
Citric acid /sodium chloride /sodium hydroxide solution	pH 4	480 mV

The potentials given refer to the standard calomel electrode at 25 °C.

4. Maintenance


Replacement of the Hydrogen cell

Unscrew the screw cap and remove the hydrogen cell.

In the recess in the electrode head you see a small O-ring. Please grease this, e. g. with silicone grease.

Take the replacement hydrogen cell and remove the protective foil. There are two openings visible. If any adhesive residue remains, remove it by rubbing with a cloth.

Insert the hydrogen cell with these two openings down into the recess.

Turn the screw cap back in.

Mark month and year on the label. Attach the sticker to the electrode head.

Place Mini-HydroFlex in water for 24 hours.

5. Troubleshooting

Potential is not displayed correctly

Cause 1: Poor contacting

Check the measuring cables on optical damages such as corrosion, cracks and sessile plugs. Replace the cable.

Check the electrode with an external reference electrode, e.g. a calomel or silver silver chloride electrode (see 3. Checking the Mini-HydroFlex).

Cause 2: No or too little hydrogen production

Did you place the new Mini-HydroFlex in water for 24 hours?

If not, wait 24 hours before you measure.

Is the duration of the hydrogen source exceeded?

If yes, please replace the hydrogen cell (see 4. Maintenance).

Cause 3: slow ion exchange, e.g. when you switch from concentrated solutions to weak concentrated solutions

Await the setting time, sometimes the balancing of concentrations lasts longer than expected. If necessary, check the potential in another electrolyte such as 1 mol/l hydrochloric acid (see 3. Checking the Mini-HydroFlex).

Cause 4: air/oxygen gets to the electrode

Avoid that gases such as air or oxygen reach the hydrogen electrode's shaft. If this happens the hydrogen is displaced or reacts off. No hydrogen potential can be set.

Fluctuating, noisy or oscillating potentials

Cause 1: Hydrogen bubbles of the hydrogen electrode itself

Bubbles are continuously leaking out of the hydrogen electrode. Some bubbles are small, others are bigger. Usually, they do not disturb your measurements.

Is a big bubble forming at the bottom of the hydrogen electrode which sticks to the vessel wall?

If possible, place the electrode further away from the vessel wall or hang the hydrogen electrode transversely into the measuring vessel.

Cause 2: Induced gases

Change the position of your gas inlet.

Induced gas bubbles which are piped by close to the hydrogen electrode can disturb the measurements and lead to fluctuating potentials.

Cause 3: Potentiostat or measuring device

Check your measuring device.

If you measure in low conductive electrolytes potentiostats and measuring devices come fast to their limits.

Deposits (e.g. red or brown) on the electrode

The hydrogen electrode is contaminated and a mixed potential is set.

Check the running time of the hydrogen source. If necessary, the hydrogen cell must be replaced as described in 4. Maintenance.

Some deposits can be wiped off with a cloth.

With 1000 grit sandpaper, you can carefully sand the bottom of the electrode.

More stubborn impurities dissolve in 1 mol / l nitric acid. The electrode must then be rinsed thoroughly in distilled water.

If that is not successful, the electrode is so badly contaminated that we suggest a new purchase.

6. Further information

For more information, please visit www.gaskatel.com.

Operating conditions

pH range:	pH -2 to pH 16
Temperature:	-20 °C to 120 °C
HF/Fluoride:	Up to 500 000 ppm
Caustic potash, caustic soda solution:	Up to 50 wt%

Technical data

Total length:	87 ± 55 mm
Shaft length:	53,5 ± 0,5 mm
Shaft diameter:	1.6 mm
Electrode shaft:	Polyether ether ketone
Electrode head:	Polypropylene
Thread:	G 1/8"
Measuring electrode:	Platinum, Palladium, PTFE
Screw cap:	Polyvinyle chloride (PVC)
Runtime hydrogen cell:	12 month

Disposal

The electrode must not be disposed of with household waste. Please dispose the electrode in accordance with local regulations in a separate collection for electrical and electronic equipment.

If sent to us free of charge (sufficiently stamped package), we will dispose of the used electrodes from our product range free of charge.

Exclusion of liability

The product is subject to continuous further development. For this reason, it is possible that parts of the instructions, technical data or images in these instructions may differ slightly from the product in front of you. The information in these instructions is only intended to clarify how the product is to be handled.

If you have any questions, please do not hesitate to contact us.

Legal claims based on these instructions cannot be asserted.