HydroFlex

Wasserstoffreferenzelektrode
Hydrogen Reference Electrode

Artikelnummer / Item number 81000 & 81010

User Manual

Version 3.2

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Lilienthalstrasse 146

D-34123 Kassel

+49 (0) 561 59190

info@gaskatel.de

www.gaskatel.com

Languages

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Operating instructions

Table of contents

Overview	2
Delivery scope	2
Declaration of conformity	2
General information	2
About these instructions	2
Intended use	2
Safety instructions	3
Use	3
Device overview	3
Operation	4
Operating conditions	10
Technical data	10
Disposal	11
Exclusion of liability	11

Overview

Delivery scope starter kit 81000

- 1 x HydroFlex hydrogen reference electrode inclusive
 - hydrogen cartridge (88010) and activation sticker (80077)
- 1 x Operating time wrench (88022)
- 1 x Test lead (89010)
- 1 x Hydrogen cartridge (88010)
- 1 x Cartridge wrench (88021)
- 1 x Activation sticker (80077)
- 1 x Manual

Delivery scope hydrogen electrode HydroFlex 81010

- 1 x HydroFlex hydrogen reference electrode inclusive
 - hydrogen cartridge (88010) and activation sticker (80077)
- 1 x Manual

Declaration of conformity

 This product complies with all EU directives applicable to this product.

General information

About these instructions

Read these instructions carefully. To ensure a long service life and reliable use, all instructions mentioned in this manual must be observed and followed. Please keep these instructions at hand.

You can access the instructions and further information at www.gaskatel.com.

Intended use

Please consider the Good Laboratory Practice (GLP).

HydroFlex is a reference electrode, exclusively intended for potential measurement during electrochemical measurements.

Safety instructions

The electrode is only intended for the above-mentioned purpose. The manufacturer is not liable for damage resulting from improper use.


Respect the following instructions, otherwise the electrode may be damaged, or the measurement results may be falsified.

When working with chemicals, all relevant safety regulations of the manufacturer and the laboratory must be respected.

Use

Device overview

Construction and dimensions

Properties of the hydrogen reference electrode HydroFlex

The hydrogen reference electrode is made entirely from plastic. Therefore, it is robust and insensitive to breakage.

It can be used in **highly concentrated alkaline solutions** as well as in **concentrated acids**, even in **hydrofluoric acid or fluoride-containing media**.

Of course, HydroFlex is **free of mercury and silver**.

Operation


1. Commissioning

Remove the protective cap from the electrode carefully.

On top of the electrode head turn the runtime clockwise to position 1 with the operating runtime key or a 3 mm hexagon key. That is corresponding to one month.

Place HydroFlex for 24 hours in distilled water. At least half of the shaft should be surrounded by water.

Mark month and year on the label. Attach the sticker to the electrode head.

After 24 hours adjust the runtime to 6 months. For correct use do not change this. The hydrogen electrode is ready for your measurements.

Please do not remove the labelling with the serial number and the sticker with the activation date.

2. Storage of hydrogen reference electrode HydroFlex

After the measurements, rinse the electrode thoroughly with water. Place the hydrogen electrode in a liquid (such as measuring solution, 1 mol/l hydrochloric acid, 1 mol/l sodium hydroxide solution, water) to keep it functional.

As soon as the running time of 6 months is reached, the hydrogen cartridge must be replaced immediately (see 4. Maintenance). Otherwise, a vacuum will form in the electrode and liquid can be soaked.

3. Checking HydroFlex

If you want to check HydroFlex, we recommend the standard calomel electrode or, if necessary, the silver silver chloride electrode.

We recommend the following solutions as measuring solution:

Hydrochlorid acid	pH 0	244 mV
Hydrochloric acid/citric acid	pH 2	362 mV
Citric acid/sodium chloride/sodium hydroxide solution	pH 4	480 mV

The potentials given refer to the standard calomel electrode at 25 °C.

4. Maintenance

Replacement of the hydrogen cell

Remove the cartridge carefully. It can be disposed of as a battery.

Take a new cartridge and grease the thread, e. g. with silicone grease.


Screw the cartridge into the head of the electrode.

You can activate the electrode now.

Complete cleaning of the hydrogen electrode


Sometimes a complete cleaning of the HydroFlex is required – especially when there are visible deposits on the electrode.

Remove the hydrogen cartridge.


Place the hydrogen electrode in 1 mol/l nitric acid for 1 hour.

Then place the entire hydrogen electrode in distilled water for 24 hours.

Dry the hydrogen electrode at about 120 °C (put the hydrogen electrode upside down) for 24 hours.

Now the hydrogen electrode is completely cleaned.

Install a new hydrogen cartridge and activate it.

5. Troubleshooting

Potential is not displayed correctly

Cause 1: poor contacting

Check the test leads on optical damages such as corrosion, cracks, and sessile plugs. Replace the test leads.

Check the electrode with an external reference electrode, e.g. a calomel or silver silver chloride electrode (see 3. Checking HydroFlex).

Cause 2: no or too little hydrogen production

Did you wait for the activation phase in water of 24 hours?

If not, wait 24 hours before you measure.

Is the duration of the cartridge exceeded?

If yes, please clean the hydrogen electrode as described above and replace the old cartridge with a new one. The new cartridge also must be activated before you can use it.

Cause 3: slow ion exchange, e.g. when you switch from concentrated solutions to weak concentrated solutions

Wait the setting time, sometimes the balancing of concentrations lasts longer than expected.

If necessary, check the potential in another electrolyte such as 1 mol/l hydrochloric acid (see 3. Checking HydroFlex).

Cause 4: air/oxygen gets to the electrode

Avoid that gases such as air or oxygen reach the hydrogen electrode's shaft. If this happens the hydrogen is displaced or reacts off. No hydrogen potential can be set.

Fluctuating, noisy, or oscillating potentials

Cause 1: hydrogen bubbles of the hydrogen electrode itself

Bubbles are continuously leaking out of the hydrogen electrode. Some bubbles are small, others are bigger. Usually, they do not disturb your measurements.

Is a big bubble forming at the bottom of the hydrogen electrode which sticks to the vessel wall? If possible, place the electrode further away from the vessel wall or hang the hydrogen electrode transversely into the measuring vessel.

Cause 2: induced gases

Change the position of your gas inlet.

Induced gas bubbles which are piped by close to the hydrogen electrode can disturb the measurements and lead to fluctuating potentials.

Cause 3: potentiostat or measuring device

Check your measuring device.

If you measure in low conductive electrolytes, potentiostats and measuring devices quickly reach their limits.

Deposits (e.g. red or brown) on the electrode

The hydrogen electrode is contaminated, and a mixed potential is set

Check the runtime of the hydrogen supply. If necessary, the hydrogen cartridge must be replaced (see 4. Maintenance).

Wipe the electrode with a cloth or clean it in 1 mol/l nitric acid and rinse it with distilled water.

If necessary, you can carry out a complete cleaning like described in 4. Maintenance.

If that is not successful, the electrode is so badly contaminated that we suggest a new purchase.

6. Further information

For more information, please visit www.gaskatel.com.

Operating conditions

pH range:	pH -2 to pH 16
Temperature:	-20 °C to 120 °C (PTFE-body only)
HF/Fluoride:	Up to 500 000 ppm
Caustic potash, caustic soda solution:	Up to 50 wt%

Technical data

Total length:	120 mm
Shaft length:	80 mm
Shaft diameter:	8 mm
Electrode shaft:	Polytetrafluorethylene (PTFE)
Electrode head:	Polytetrafluorethylene (PTFE)
Measuring electrode:	Platinum/Palladium/PTFE
Hydrogen cartridge:	Polycarbonate (PC)
Runtime hydrogen cell:	Recommended 6 months

Disposal

The electrode must not be disposed of with household waste. Please dispose of the electrode in accordance with local regulations in a separate collection for electrical and electronic equipment.

If sent to us free of charge (sufficiently stamped package), we will dispose of the used electrodes from our product range free of charge.

Exclusion of liability

The product is subject to continuous further development. For this reason, it is possible that parts of the instructions, technical data or images in these instructions may differ slightly from the product in front of you. The information in these instructions is only intended to clarify how the product is to be handled.

If you have any questions, please do not hesitate to contact us.

Legal claims based on these instructions cannot be asserted.