HydroFlex

Wasserstoffreferenzelektrode

Artikelnr.: 81010 , 81000

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
Operating Instructions

Intended Use

This hydrogen electrode is a reference electrode intended only for potential measurement during electrochemical measurements.
Note: The hydrogen source of this reference electrode must be activated before first use.

Safety Information

The electrode is intended only for the purpose above. The manufacturer is not liable for damage resulting from improper use.
Observe the following instructions, otherwise the electrode may be damaged or measurements results be falsified.
When working with chemicals, all relevant safety regulations of the manufacturer and the laboratory must be observed.

Commissioning

Provide a solidly standing beaker with distilled or deionized water.

Remove the protective cap from the electrode carefully.


On top of the electrode head turn the runtime clockwise to position 1 with a 3 mm hex screw driver. That is corresponding to one month.

Place Hydroflex for 24 hours in the water. At least half of the shaft should be surrounded by water.

Important: Mark the activation month and year on the enclosed sticker and attach it to the HydroFlex.

After 24 hours adjust the runtime to 6 months. For correct use do not change this. The Hydrogen electrode is ready for your measurements.

Storage

After the measurements, rinse the electrode thoroughly with water. Place the hydrogen electrode in a liquid (such as measuring solution, 1 mol / l hydrochloric acid, 1 mol / l sodium hydroxide solution, water) to keep it functional.

Checking the HydroFlex

If you want to check HydroFlex, we recommend the standard calomel electrode or, if necessary, the silver silver chloride electrode. We recommend the following solutions as measuring solution:

Hydrochlorid acid	pH 0	244 mV
Hydrochloric acid/citric acid	pH 2	362 mV
Citric acid/sodium chloride/sodium hydroxide solution	pH 4	480 mV

The potentials given refer to the standard calomel electrode at 25°C.

Runtime of Hydrogen Cartridge

Once the end of runtime has been reached, HydroFlex must no longer be stored in liquid. Replace the cartridge and activate the reference electrode as described in commissioning.

Remove the cartridge.

Take the new cartridge and grease the thread.

Screw the cartridge back into the head of the electrode.

You can activate the electrode now.

Troubleshooting

Make sure all cables are connected correctly. Check if the connections are corroded or damp. Note the lifetime of the hydrogen source! If necessary, replace the cartridge. Atmospheric oxygen distorts the hydrogen potential. Therefore leave HydroFlex in air for a maximum of 5 seconds.

Further Informations

More information can be found on our website www.gaskatel.com.

Disposal

The electrode must not be disposed of as household waste. Please dispose of the electrode in a separate collection for electrical and electronic equipment according to local regulations. This also applies to the hydrogen cells. If you have further questions, please contact us.