Mini-HydroFlex

Wasserstoffreferenzelektrode

Artikelnr.: 81020

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
Operating Instructions

Intended Use

The mini-HydroFlex is a reference electrode intended only for potential measurement during electrochemical measurements.

Safety Information

The electrode is intended only for the purpose above. The manufacturer is not liable for damage resulting from improper use. Observe the following instructions, otherwise the electrode may be damaged or measurement results may be falsified. When working with chemicals, all relevant safety regulations of the manufacturer and the laboratory must be observed.

Commissioning / Replacement of Hydrogen cells

Unscrew the black cap and remove the hydrogen cell.

In the recess in the electrode head you see a small o-ring. Please grease this, e.g. with silicone grease.

Take the replacement hydrogen cell and remove the protective foil. There are two openings visible. If any adhesive residue remains, remove it by rubbing with a cloth.

Insert the hydrogen cell with these two openings down into the recess.

Turn the cap back in.

Mark month and year on the label. Attach the sticker to the electrode head.

Place Mini-HydroFlex in water for 24 hours.

Please do not remove the labelling with the serial number and the sticker with the activation date.

Storage

After the measurements, rinse the electrode thoroughly with water. Place the hydrogen electrode in a liquid (such as measuring solution, 1 mol / l hydrochloric acid, 1 mol / l sodium hydroxide solution, water) to keep it functional.

Runtime of Hydrogen source

Pay attention to the runtime of the hydrogen source. This is **12 month**. Once the end of runtime has been reached, Mini HydroFlex must not longer be stored in liquid. Replace the hydrogen source immediately as described in commissioning / Replacement of hydrogen cell.

Checking the HydroFlex

If you want to check HydroFlex, we recommend the standard calomel electrode or, if necessary, the silver silver chloride electrode. We recommend the following solutions as measuring solution:

Hydrochlorid acid	pH 0	244 mV
Hydrochloric acid/citric acid	pH 2	362 mV
Citric acid/sodium chloride/sodium hydroxide solution	pH 4	480 mV

The potentials given refer to the standard calomel electrode at 25°C.

Troubleshooting

Make sure all cables are connected correctly. Check if the connections are corroded or damp. Note the lifetime of the hydrogen source! If necessary, replace them. Atmospheric oxygen distorts the hydrogen potential. Therefore leave Mini-HydroFlex in air for a maximum of 5 seconds.

Further Informations

More information can be found on our website www.gaskatel.com.

Disposal

The electrode must not be disposed of as household waste. Please dispose of the electrode in a separate collection for electrical and electronic equipment according to local regulations. This also applies to the hydrogen cells. If you have further questions, please contact us.