pHydrunio

Wasserstoff-pH-Elektrode

Artikelnr.: 84100

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
Operating Instructions

Intended Use

The Hydrogen-pH-electrode pHydrunio intended only for measurements of pH in a range from pH- 2 to 16. pHydrunio is delivered with activated hydrogen cells. But you must fill in the inner buffer solution (commissioning).

Safety Information

The electrode is intended only for the above purpose. The manufacturer is not liable for damage resulting from improper use. Observe the following instructions, otherwise the electrode may be damaged or measurements results be falsified. When working with chemicals, all relevant safety regulations of the manufacturer and the laboratory must be observed.

Commissioning

Provide a solidly standing beaker with calibration buffers pH 7.

Remove the red cap on the electrode head. 2 openings are visible. Fill 5 ml of the inner buffer pH 7 into the larger one with syringe and cannula. Screw the red cap back on.

Place pHydrunio in calibration buffer pH 7 for 24 hours.

Please do not remove the labelling with the serial number and the sticker with the activation date.

Storage

After the measurements rinse pHydrunio thoroughly with water. Store pHydrunio in calibration buffer pH 7 if you do not measure.

Calibration

First, the pH-electrode is checked in the calibration buffer pH 7 (Gaskatel). Note the temperature as well as the measured voltage. Then rinse the pH-electrode and the temperature probe well with water and place it in the calibration buffer with pH 2.0. Note the temperature as well as the measured voltage.

Calculation measured potential to – pH-value

The measured voltage U of can be converted into the pH value as follows:

$$\text{pH} = \frac{-E \text{ in Volt}}{273,15 + \text{Temperatur in } ^\circ\text{Celsius}} \cdot 5039,1 + 7$$

Runtime of Hydrogen cells


Pay attention to the runtime of the hydrogen source. This is 12 month. Once the end of runtime has been reached, pHydrunio must no longer be stored in liquid. Replace the hydrogen cells immediately.

There are two black caps on the sides of the electrode head. Please unscrew them, too. Remove the hydrogen cells.

In the recess in the electrode head you see a small o-ring. Please grease this, e.g. with silicone grease.

Take the replacement hydrogen cell and remove the protective foil. There are two openings visible. If any adhesive residue remains, remove it by rubbing with a cloth.

Insert the hydrogen cell with these two openings down into the recess..

Turn the cap back in.

Mark month and year on the label. Attach the sticker to the electrode head.

Place the hydrogen pH electrode in calibration buffer pH 7 for 24 hours .

Troubleshooting

Make sure all cables are connected correctly. Check whether the connections are corroded or damp. Note the runtimes of the hydrogen button cells! If necessary, replace them. Check whether you are measuring in the correct calibration liquid. Replace the reference electrolyte. Atmospheric oxygen distorts the hydrogen potential. Therefore leave pHydrunio in air for a maximum of 5 seconds.

Further informations

More information can be found on our website www.gaskatel.de.

Disposal

The electrode must not be disposed of as household waste. Please dispose of the electrode in a separate collection for electrical and electronic equipment according to local regulations. If you have further questions, please contact us.